

# Ultrasonic Cleaner ULTRA NEUTRAL

## Description

The liquid ultrasonic cleaner ULTRA NEUTRAL is a concentrate und universally suitable for gold, silver, steel and non-ferrous metals. ULTRA NEUTRAL meets the legal requirements with over 80% biodegradability according to OECD methodology. If properly used, it does not cause any problems with wastewater treatment according to past experiences.

## Storage

Store it closed. No further conditions.

## Risks/Disposal

Used solutions should be disposed of by a specialist company who will treat it before discharging. Check with the local water authority for permitted discharge levels.

## Operating conditions

Preparation	3 – 5%
Temperature	60 – 80°C
pH value	7 - 8

**Please note our  
Safety Data Sheet!**

## Form of delivery

- Concentrate, liquid

## Bath preparation

Per liter ready to use bath needed: 30– 50 ml  
Ultrasonic cleaning concentrate ULTRA  
NEUTRAL.

## Special hints

With new application, the heated up cleaning liquid must be vented through examine “without goods” approx. 1-5 minutes.

The necessary cleaning time varies depending on the type of goods and the degree of soiling. The times are to be determined by tests, cleaning times over 3 minutes are pointless.